



CPA/CPB/CPC Series

- Resolution 0.01g, 0.1g



CPA102-CPA2002



CPB122-CPB3002



CPC1002-CPD10001

| Features                        | Benefits                                                 |
|---------------------------------|----------------------------------------------------------|
| Economical and reliable         | Helps you save money with confidence                     |
| Large stainless steel pan       | Allows for easy weighing of large objects                |
| Leveling feet with level bubble | Makes it level even on uneven working surfaces           |
| Below balance weighing          | Determines density of large pieces                       |
| AC-DC available                 | Able to operate with rechargeable battery (not included) |
| LCD display with backlight      | Allows for easy reading                                  |
| Multiple weighing units         | Allows for flexibility of weighing                       |
| Double TARE key                 | Easy for work to operate                                 |

| Model             | CPA<br>302                                               | CPA<br>622 | CPA<br>1202 | CPB<br>602             | CPB<br>1002 | CPB<br>2002 | CPB<br>3002 | CPC<br>2002            | CPC<br>3002 | CPC<br>5002 | CPD<br>6001 | CPD<br>10001 |
|-------------------|----------------------------------------------------------|------------|-------------|------------------------|-------------|-------------|-------------|------------------------|-------------|-------------|-------------|--------------|
| Draft Shield      | With Plastic Draft Shield                                |            |             |                        |             |             |             | No Draft Shield        |             |             |             |              |
| Capacity          | 300g                                                     | 620g       | 1200g       | 600g                   | 1000g       | 2000g       | 3000g       | 2000g                  | 3000g       | 5000g       | 6000g       | 10kg         |
| Resolution        | 0.01 g                                                   |            |             |                        |             |             |             |                        |             |             | 0.1 g       |              |
| Repeatability     | ± 0.01 g                                                 |            |             |                        |             |             |             |                        |             |             | ± 0.1 g     |              |
| Linearity Error   | ± 0.02 g                                                 |            |             |                        |             |             |             |                        |             |             | ± 0.1 g     |              |
| Response Time     | ≤ 2 sec                                                  |            |             |                        |             |             |             |                        |             |             |             |              |
| Calibration       | Calibration with External Weight                         |            |             |                        |             |             |             |                        |             |             |             |              |
| Weighing Units    | g, ct, lb, oz, ozt, dwt, kg, tmR, GSM, PKT etc(17 units) |            |             |                        |             |             |             |                        |             |             |             |              |
| Construction      | Cast aluminum, static resistant ABS, stainless steel pan |            |             |                        |             |             |             |                        |             |             |             |              |
| Display           | Large, easy to read LCD display with backlight           |            |             |                        |             |             |             |                        |             |             |             |              |
| Pan Size          | Round Pan (115mm)                                        |            |             | Square Pan,157 x 136mm |             |             |             | Square Pan,192 x 171mm |             |             |             |              |
| Body dim.         | 304x 204 x 300 mm                                        |            |             |                        |             |             |             | 304x 204 x 90 mm       |             |             |             |              |
| Battery           | Rechargeable battery                                     |            |             |                        |             |             |             |                        |             |             |             |              |
| Power             | 110v - 240v, 50HZ / 60HZ                                 |            |             |                        |             |             |             |                        |             |             |             |              |
| Operational Temp. | 10 °C to 40 °C                                           |            |             |                        |             |             |             |                        |             |             |             |              |
| Connectivity      | Configurable RS 232 interface                            |            |             |                        |             |             |             |                        |             |             |             |              |