

OVERHEAD STIRRER

Direct drive overhead Stirrer with speed regulator

Overhead stirrer is used to mix large volume of liquid and gel in laboratory. Our laboratory overhead stirrer comes with heavy duty motor and drill chuck.

Standard fittings are heavy cast iron base, height adjustable stand and speed regulator box. Standard model comes with 1/20 HP motor; in addition, 60-minute timer can also be fitted as an option with added cost.



TECHNICAL SPECIFICATIONS

Model	STXOHS12
RPM	4000
Motor	1/20 HP
Speed regulator	Speed regulator in separate box
Standard fittings	• Height adjustable Stand • Cast iron base • Rotator
Power supply	220 Volts 50Hz
Optional	• 0 - 60 minutes timer • Digital RPM display