

STERICOX

SEDIMENTATION SHAKER

Sedimentation shaker for 5 or 6 cylinders

STERICOX™ sedimentation shaker is designed for smooth analysis of sedimentation value of samples under 50 or 100 ml cylinders. These shakers are made for 5 or 6 cylinders; cylinder carrier (holding tray) is made either wooden or acrylic. These units feature smooth operation and movement at an angle of 30° from horizontal. For shaking / oscillation a speed regulator is fitted at front side of the cabinet; by rotating it, user can set the RPM as required. A digital timer with buzzer is an optional accessory.



It can be fitted on extra price which alerts user when cycle (sedimentation period) ends. It is light weight in design, cost effective and offers long years of trouble-free working performance.

TECHNICAL SPECIFICATIONS

Model	STXVDR13	STXVDR13D
Platform size	13 x 13 inch with spring holder	13 x 13 inch with spring holder
RPM	150 to 200 rpm	150 to 200 rpm
RPM display	NO	LED Display
Timer	Mechanical Timer 0-60 minutes	Digital Timer 0-60 minutes
Power supply	220 Volts 50 Hz	220 Volts 50 Hz

STERICOX

STERICOX INDIA PRIVATE LIMITED

Address : Ground Floor A-106, KH No. 186
Street No. 15. South Gamri. Delhi 110053. India
Call Us : +918130134164
Email Us : info@stericox.com
Web : www.stericox.com